## IFU

(Instruction for use)

### Manípulo do Análogo Digital / Digital Analogue Handle / Manípulo del Análogo Digital

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to incorporate clinical, mechanical and biological concepts, so for best results it should only be used with original associated components in accordance with Emfils instructions and recommendations. All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthetists and Oral and Maxillofacial Surgeons.

**DESCRIPTION**

The handles make it easy to handle, transport and fit the Emfils digital analog into the prototyped model. The raw material used in the EMFILS assembler is ASTM F899 stainless surgical steel.

The handles are presented according to the table below:


Sales code	Conection Type				Thread	
	CM/LZ	HE	HI	Transmucosal		
				Multi-Unit		Mini 4,8
MAD-1414					X	M1,4
MAD-1616	X					M1,6
MAD-2020		X	X	X		M2,0

**APPLICATION**

The handles are indicated for transport and installation of the digital analog in the cavity of the prototyped working model.

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It is contraindicated for use on side loads. Apply load only along key axis.

**HANDLING**

Screw in the respective digital analogue and bring the assembly to the cavity of the prototyped model. Always press along the long axis. To remove the digital analogue pull the assembly, always keeping the shaft.

All procedures must be performed by qualified professionals and the handling of the product is inherent to the training of the professional, who will be responsible for using the instruments in the best possible way, following the instructions of Emfils for each instrument.

Important points in use:

- Avoid side loading.
- The product has no torque control.
- Use of excessive force may damage the product.

**PRESENTATION**

The product is individually packed in tyvek® type surgical paper packaging and transparent film display. The instrument undergoes a decontamination process that eliminates any dirt, but the cleaning process described above should be followed at first use.

**HYGIENIZATION AND STERILIZATION**

The handle is reusable, supplied non-sterile and must be sanitized and sterilized before and after use when used in the surgical procedure, following the recommendations of the document available at the link below:

[Novo Colosso's Surgical Kit Cleaning Manual](#)

**PRECAUTIONS**

- Do not use if packaging is broken;
- Avoid side loading on the product;
- The product has no torque control.
- Before each procedure check the condition of the instruments, always respecting their useful life. Replace the instruments in case of damage, erased markings, deformation and wear;
- Excessive force may compromise efficiency and damage the product or component connected to it;
- Always use the product sequence Emfils. The use of instruments and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Emfils implant system and disclaims any product warranty.

**ADVERSE EFFECTS**

Adverse effects will only occur if the choice of instrument is inadequate, which may cause damage to the component to be installed.

**STORAGE CONDITIONS**

The product should be stored in a dry place and in its original packaging at room temperature and not exposed to direct sunlight. Improper storage may influence the original packaging and product characteristics.

**DISPOSAL OF MATERIALS**

The professional should always evaluate the condition of the instruments after each use, and if detected that the instrument presents damage, erased markings, deformation and wear, should be discarded in accordance with current regulations, in compliance with the environmental requirements and the criteria of potentially hazardous products, taking into account the different levels of contamination, physical hazards and selective collection

**PRAZO DE VALIDADE**

Indicated on product label

Not all products are available in all countries.

Please contact your authorized distributor.

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