## IFU

(Instruction for use)

### Expansor com Rosca Guiada / Guide Expander with Thread / Expansor con Rosca Guiada

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

**DESCRIPTION**

The bone expander screw Guide form a set of four parts, available in short and long lengths, and is indicated for the opening of the surgical cavity through the bone growth process and lateral compression to the guided surgery technique before installation implants Ø3,3mm / 3,5mm / Ø4,0mm / Ø4,1mm with thin bone bed type III and IV, according to the bone classification of lekholm and zarb (1985). Bone classification of Lekholm and zarb (1985):

bone type I = homogeneous compact bone.

bone type II = compact bone with trabecular bone nucleus.

bone type III = thin layer of cortical bone around the dense trabecular bone of good resistance.

bone type IV = thin layer of cortical bone around the trabecular bone with low density.

At its active end, the product contains a thread that facilitates and guides the insertion of the part into the bone. At the other end the threaded bone expander has a universal square fit of 3.9 mm, which connects to the torque wrench gear. The laser markings take into account the actual implant length, with an additional 11 mm length to compensate for the soft tissue thickness and height of the washer and the Drill Guide.

The Bone Expanders with Screw Guide are presented in the table below:

Guide Bone Expanders with Thread			
Sales Code	Diameter	Model	Total Length (mm)
ECG-3200	3,2	Short	27
ECG-3800	3,8	Short	27
ELG-3200	3,2	Long	33
ELG-3800	3,8	Long	33

Bone Guide Expanders receive the application of a thin biocompatible carbon coating that imparts dark coloration to the cutting surface. This film has as main purposes:

- Increased resistance to oxidation;
- Reduced friction between drill and bone;
- Increased wear resistance;
- Reduction of the heat generated in the bone during the osteotomy;
- Optimization of the visualization of height markings;

The raw material used is stainless steel tempered steel, according to ASTM F899.

**APPLICATION**

The bone with expanders Guide thread guides are surgical instruments used to perform bone expansion technique that prepares the surgical alveolus for installing implants in cases with little rim thickness. They are only indicated where the thickness of the alveolar crest has a minimum of 3 mm.

**CONTRAINDICATIONS**

Use in dense bone beds type I and II (according to the bone classification of lekholm and zarb (1985)) is contraindicated.

**HANDLING**

Once sterilized, the Bone Screw Expander Guide should be handled only in a sterile environment, by properly trained professionals, and in appropriate attire at the time of implant surgery.

All procedures must be performed by qualified professionals, and the handling of the product is inherent in the training of the professional, who will be responsible for using the instruments in the best possible way, respecting the indications of the Emphasis for each instrument.


To conform the surgical alveolus with Guide Bone Bone Expander, after fixing the Surgical Guide into the oral cavity, select the Emfils Guide Drill Guide corresponding to the expander that will be used and according to the washer model. Each Expansion Guide has a specific Emfils Guide Drill Guide. Snap the end of the Drill Guide into the washer of the Surgical Guide, holding the cable firmly. Insert the Expansion Guide into the Torque Wrench, after making insertion movements, clockwise by lightly pressing the Torquimetro Ratchet Wrench, with a maximum torque of 60 N.cm. Insert the Expander Guide until the marking for the chosen implant or until the end of the stop according to virtual planning. No irrigation is required.

Comments:

1. Always keep the instrument in the direction of insertion.
2. The pressure on the Bone Expander with Thread Guide should be controlled.

**SANITATION AND STERILIZATION**

Bone Screwdrivers Guide are reusable, supplied non-sterile and must be sanitized and sterilized before and after use, following the recommendations in the document available at the link below:

**PRESENTATION**

The Bone Screwdrivers Guide Guide are individually packaged in tyvek® surgical paper packaging and clear film display. The instruments are also sold in kits mounted on the case in an identified position, and are affixed by silicone rings or brackets. The kit is guaranteed by exclusive seals in both the length and width of the case. The instruments undergo a decontamination process that eliminates any type of dirt, but the hygiene process, described above, must be followed as soon as it is first used.

**PRECAUTIONS**

- Do not use if the packaging has been tampered with;
- This product should be used sterile. Do not sterilize under dry heat (greenhouse), always use steam autoclave and pressure;
- Chemical sterilization is not recommended;
- Do not autoclave the product in its original packaging;
- When using Bone Screw Extenders Guide should be careful not to cause fractures in the lateral bone walls of the surgical alveolus;
- Care must be taken in the hygienization because it is a sharp metal instrument, so as not to suffer any type of physical trauma;
- Make sure to use the diameter compatible with the procedure and with the diameter of the planned implant;
- Instruments should not be re-established;
- Due to the reduced opening of the mouth in the posterior region, the use of the instrument may be difficult;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient).
- Make sure you have all the necessary instruments to perform the surgery according to surgical planning;
- Before each procedure, check the conditions of the instruments in the first, always respecting their useful life. Replace instruments in case of damage, missing markings, deformation and wear;
- The use of excessive force may compromise the procedure;
- Always use the product sequence in file. The use of instrumentation and / or prosthetic components from other manufacturers does not guarantee the perfect function of the Implant System and exempts any warranty of the product.

#### ADVERSE EFFECTS

Adverse effects will only occur if the choice of instruments is inadequate, which can lead to damage to the prosthetic sequence of the installed implant.

#### STORAGE CONDITIONS

The product should be stored in a dry place in its original packaging at room temperature and not exposed to direct sunlight. Incorrect storage may influence the original characteristics of the carton and the product.

#### DISPOSAL OF MATERIAL

The operator must always evaluate the condition of the instrumentation after each use and it must be discarded in accordance with current regulations, complying with the environmental requirements and the criterion of potentially dangerous products, and if the instrument is damaged, markings erased, deformed and worn, taking into account the different levels of contamination, physical risks and selective collection.

#### LIFE CYCLE

Bone Screwdrivers Guide should be discarded when they lose their functionality or the user detects excessive use.

#### SHELF LIFE

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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