IFU


(Instruction for use)

Protetor / Protective / Protector

PORTUGUÊS

ENGLISH

ESPAÑOL

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This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The protectors are part of the the prosthetic components family for transient use that are installed on the base abutment to protect the entrance of the thread of the same against particles that can have access to this region during the preparation of the prosthesis and / or cicatrization of the peri-implant tissues, restricting the prosthesis screw access. The protector has at one tip, a fitting for the hexagonal 1.2 mm driver, that allows its installation. At the other tip, the protector has a compatible geometry and thread for installation on the corresponding base abutment.

The raw material used is the titanium alloy (grade 5), according to ASTM F136.

The protectors are presented in the table below:

Sales Code	Base Abutment Application				Base Abut. Diameter	Thread	Max. Torque
	Multi	Mini 4,0-5,0	Mini 4,8	Angled			
PPT-4004	X				4,0	M2	15 N.cm
PPT-5004	X				5,0	M2	15 N.cm
PM-4004		X			4,0	M1.6	15 N.cm
PM-5004		X			5,0	M1.6	15 N.cm
PMIP-4803			X	X	4,8	M1.4	15 N.cm

APPLICATION

The protectors are indicated for short period use to protect the entry of the base abutment thread during the final confection of the prosthesis and / or healing of the peri-implant tissues. They must be installed in the mouth on the corresponding base abutment, so always check the mechanical compatibility between the components. Make sure that the guard has been tightened to the required torque according to the table above.

CONTRAINDICATIONS

This product is not for customization.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.emfils.com.br/ifu>

HANDLING

Collect the protector using the length of 1.2 mm hexagonal driver that best suited to the work you wish to perform. At the time of installation, be sure to align the screw along the corresponding base abutment shaft. The component must be perfectly seated to the base abutment after application of the necessary torque. The maximum torque value to be applied to the protector must be respected.

PRESENTATION

Any protector is indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Any type of protective is for single use only.
- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause

damage to the material and make the system unusable;

- It is recommended that the protective is remove from the patient’s mouth within 6 months after installation;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the protector along the component shaft, thus preventing thread damage or system locking;
- The component must be perfectly seated to the base abutment after the necessary torque has been applied to the protective. This check is recommended using periapical radiographs with the parallelism technique, when applicable.
- Instruct the patient that if he or she notices that the protective will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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