IFU


(Instruction for use)

Pilar Calcinável CoCr / CoCr Burnout Abutment / Pilar Calcinable CoCr

PORTUGUÊS

ENGLISH

ESPAÑOL

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Código de Venda	Tipo Conexão								Tipo Prótese		Diâmetro Implante ou Transm.	Altura Cinta	Altura Ci- mentação	Parafuso
	CM	LZ	HE	HI	Transmucoso				Unit.	Múlt.				
					Multi	Mini 4,0-5,0	Mini 4,8	Angul.						
CLNE-3511			X						X		3,5	---	9,0	M2
CLNE-4011			X						X		4,0	---	9,0	M2
CLNE-5011			X						X		5,0	---	9,0	M2
CLNI-3511		X							X		3,5	---	9,0	M1.6
CLNI-4011		X							X		4,0	---	9,0	M1.6
CLNI-5011		X							X		5,0	---	9,0	M1.6
CLNNR-4010					X				X		3,5-4,0-5,0	---	10,0	M2
CLNNR-5010					X				X		3,5-4,0-5,0	---	10,0	M2
CLRE-3511			X							X	3,5	---	9,0	M2
CLRE-4011			X							X	4,0	---	9,0	M2
CLRE-5011			X							X	5,0	---	9,0	M2
CLRI-3511		X								X	3,5	---	9,0	M1.6
CLRI-4011		X								X	4,0	---	9,0	M1.6
CLRI-5011		X								X	5,0	---	9,0	M1.6
CLRMP-4010						X				X	3,5-4,0-5,0	---	10,0	M1.6
CLRMP-5010						X				X	3,5-4,0-5,0	---	10,0	M1.6
CLRNM-3507	X								X		3,5-4,0-5,0	1,0	9,0	M1.6
CLRNM-3515	X								X		3,5-4,0-5,0	1,5	9,0	M1.6
CLRNM-3530	X								X		3,5-4,0-5,0	3,0	9,0	M1.6
CLRNM-3545	X								X		3,5-4,0-5,0	4,5	9,0	M1.6
CLRNM-3560	X								X		3,5-4,0-5,0	6,0	9,0	M1.6
CLRNR-4010					X					X	3,5-4,0-5,0	---	10,0	M2
CLRNR-5010					X					X	3,5-4,0-5,0	---	10,0	M2
CPPLE-3315				X						X	3,3-3,5	---	15,0	M2
CPPLE-4015				X						X	4,0	---	15,0	M2
CPPLE-5015				X						X	5,0	---	15,0	M2
CPPLEH-3315				X					X		3,3-3,5	---	15,0	M2
CPPLEH-4015				X					X		4,0	---	15,0	M2
CPPLEH-5015				X					X		5,0	---	15,0	M2
LMPP-4810							X	X		X	3,5-4,0-5,0	---	10,0	M1.4

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfil instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The burnout abutments or cylinders with Cobalt-Chromium alloy for micro-casting are

part of the prosthetic components family and are used as an intermediate between the implant or the base abutment and the prosthesis (crown). The product is available with different fittings, depending on the type of implant or base abutment used that can have a rotational or non-rotational characteristic.

It is made of a high melt Cobalt - Chromium base alloy and receives a cylindrical cover in burnout polymer for use in the laboratory in the phases of waxing and the micro-casting of screwed prostheses infrastructures on the implant or base abutment. It has a cylindrical shape that can be prepared in the best way

to meet the need of the work, as long as it also meets the instructions described without compromising the product.

It allows micro-casting of noble, semi-noble and alternative alloys with a melting point of up to 1250° C (2282° F). Density 8.4 g / cm³ (520 lb / ft³). Coefficient of thermal expansion 13 µm / m-K.

The burnout abutments or cylinders in CoCr are presented in the table below:

Sales Code	Connection Type								Prosthesis Type		Implant Diameter or Base Abut.	Collar Height	Cimentation Height	Screw
	MT (CM)	ZL (LZ)	EH (HE)	IH (HI)	Base Abutment				Unit.	Múlt.				
					Multi	Mini 4,0-5,0	Mini 4,8	Angled						
CLNE-3511			X						X		3,5	---	9,0	M2
CLNE-4011			X						X		4,0	---	9,0	M2
CLNE-5011			X						X		5,0	---	9,0	M2
CLNI-3511		X							X		3,5	---	9,0	M1.6
CLNI-4011		X							X		4,0	---	9,0	M1.6
CLNI-5011		X							X		5,0	---	9,0	M1.6
CLNNR-4010					X				X		3,5-4,0-5,0	---	10,0	M2
CLNNR-5010					X				X		3,5-4,0-5,0	---	10,0	M2
CLRE-3511			X							X	3,5	---	9,0	M2
CLRE-4011			X							X	4,0	---	9,0	M2
CLRE-5011			X							X	5,0	---	9,0	M2
CLRI-3511		X								X	3,5	---	9,0	M1.6
CLRI-4011		X								X	4,0	---	9,0	M1.6
CLRI-5011		X								X	5,0	---	9,0	M1.6
CLRMP-4010						X				X	3,5-4,0-5,0	---	10,0	M1.6
CLRMP-5010						X				X	3,5-4,0-5,0	---	10,0	M1.6
CLRNM-3507	X								X		3,5-4,0-5,0	1,0	9,0	M1.6
CLRNM-3515	X								X		3,5-4,0-5,0	1,5	9,0	M1.6
CLRNM-3530	X								X		3,5-4,0-5,0	3,0	9,0	M1.6
CLRNM-3545	X								X		3,5-4,0-5,0	4,5	9,0	M1.6
CLRNM-3560	X								X		3,5-4,0-5,0	6,0	9,0	M1.6
CLRRN-4010					X					X	3,5-4,0-5,0	---	10,0	M2
CLRRN-5010					X					X	3,5-4,0-5,0	---	10,0	M2
CPPLE-3315				X						X	3,3-3,5	---	15,0	M2
CPPLE-4015				X						X	4,0	---	15,0	M2
CPPLE-5015				X						X	5,0	---	15,0	M2
CPPLEH-3315				X					X		3,3-3,5	---	15,0	M2
CPPLEH-4015				X					X		4,0	---	15,0	M2
CPPLEH-5015				X					X		5,0	---	15,0	M2
LMPP-4810							X	X		X	3,5-4,0-5,0	---	10,0	M1.4

APPLICATION

The burnout abutments or cylinders with a CoCr termination are suitable for unitary (non-rotational fitting) or multiple (rotational fitting) prostheses, made from the waxing and micro-casting of screwed prosthesis infrastructures on the base abutment or implant, according to available interocclusal space, existing gingival height and three-dimensional position of the implant. Indicated to be used in immediate rehabilitation procedures, in the maxilla or

mandible. Before finishing the clinical work, make sure that the screw has been tightened with the necessary torque.

CONTRAINDICATIONS

Contraindicated in cases where the interocclusal space is insufficient and / or the three-dimensional position of the implant is unsatisfactory.

MAGNETIC RESONANCE IMAGING (MR) – SAFETY INFORMATION

These products are made from a metallic material that can be affected by MR energy.

For more information, see “Safety Information | Magnetic Resonance Imaging (MR)” at <https://www.emfil.com.br/ifu>

HANDLING

After installing the implant or base abutment in the mouth, transfer its position through molding with the help of the corresponding transfer, according to appropriate techniques. Screw the burnout or cylinder over the analogue and customize according to the available interocclusal space. Carry out the waxing laboratory process, which must be at least 0.5 mm thick and can be reduced to 0.3 mm after the micro-casting. The obtained pattern should be included in an alloy compatible coating that will be used in the micro-casting process. The micro-casting alloy chosen should be an alloy also based on Cobalt - Chromium and compatible with the metal material of the abutment termination and with the aesthetic material that is planned for the prosthesis. The original shape must be preserved without the presence of positive metal bubbles. The metal structure is thus obtained which will serve as the basis for the acrylization or porcelain application. Checks and tests of passivity and adaptation of the prosthesis structure must be performed. During the finishing and polishing procedures, the prosthetic interface of the component must be insulated by the respective polishing protector. During the alloy finishing process, do not use corrosive materials, such as those containing iron particles. It is also recommended to machine the alloy in low rotation and under-refrigeration.

It is recommended to use stainless steel screw for lab work and the titanium alloy screw for fixing in the mouth, always taking into account the recommended torque.

PRESENTATION

The CoCr Burnout Abutment are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- **WARNING:** For micro-casting, always use Cobalt - Chrome alloy with lower melting temperature than above. Other metal alloys do not ensure proper bonding and may further deform the prosthetic interface of the component;
- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Any type of abutment is for single use only and should not be re-sterilized;
- Reusing this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- This product should not be used as a definitive plastic element. The micro-casting process is imperative before any installation in the mouth;
- The choice of the micro-casting alloy in the laboratory should be made by the dental surgeon in order to minimize the possible adverse effects for each patient;
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Special care should be taken during waxing for what the interocclusal space is not insufficient or excessive;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;
- The component must be perfectly seated to the implant, base abutment or analog after the necessary torque has been applied to

the screw. This check is recommended using periapical radiographs with the parallelism technique, when applicable.

- When using a unitary restoration, it is medically advisable to maintain a minimum foundation height of 4 mm above the gum;
- It is recommended that the final prosthesis be hygienical before installation in the mouth;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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