## IFU


(Instruction for use)

Base de Titânio Ti-Base 4 CEREC® / Titanium Base Ti-Base 4 CEREC® / Base de Titanio Ti-Base 4 CEREC®

PORTUGUÊS

ENGLISH

ESPAÑOL

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This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results, it should be used only with original associated components According to emfils instructions and recommendations. All procedures must be performed by qualified professionals with specific training in CAD/CAM technology, in the Sirona Dental CAD/CAM System for digitizing and in the production of custom prostheses.

**DESCRIPTION**

**TITANIUM BASE TIBASE 4 CEREC® EMFILS:** The Titanium Base TiBase 4 CEREC® Emfils is an intermediate prosthetic abutment Indicated for unitary cementation or bolted restorations on dental implants installed in the jaw. It is manufactured in Titanium alloy, according to ASTM F136, to be installed between the implant and the prosthesis, which is manufactured from the Sirona Dental CAD/CAM System for digitizing and production of custom prostheses. The Titanium Base TiBase 4 CEREC® Emfils Has a cylindrical shape with indexator guide (Notch) to fit the ScanBody, prosthesis or Coping. It has a cemable area height of 4.7 mm and with variations in height of Transgingival (see table). The geometry of the cemable area of the titanium Base TiBase 4 CEREC® Emfils is compatible with the blocks Incoris, the Zi Meso Or Disilicate Available for CEREC®. The Titanium Base TiBase 4 CEREC® Emfils is available for Use In implants Emfils of Interfaces Morse Taper, External Hexagon and Internal Hexagon. The fixing screws do not accompany the titanium bases TiBase 4 CEREC® Emfils, are sold separately. Note: Single use Only.

**SCANPOST 4 CEREC® EMFILS:** Scanpost is an extender device about implants For scanning, with plug equivalent to the TiBase 4 CEREC® Emfils that can be used Intra-oral Or Extraorally (Benchtop scanner). It is manufactured in titanium alloy, according to ASTM F136 Standard, has a cylindrical shape with indexator guide (Notch) to fit the ScanBody And is bolted In Implant to the prosthetic connection. The ScanPost 4 CEREC® Emfils can be used up to 50 times. If the screws are lost or damaged, they must be replaced. The ScanPost 4 CEREC® Emfils is available for use in emfils implants of Morse Taper, External Hexagon and Internal Hexagon Interfaces. The fixing screws do not accompany the ScanPost 4 CEREC® Emfils, are sold separately.

⚠ Note: the Scanpost 4 CEREC® Emfils should not be used as a definitive prosthetic abutment!

⚠ Note: For greater scanning accuracy, Emfils recommends the acquisition of Original Scanbodies Dentsply Sirona.

The abutments 4 CEREC® Emfils are presented in the table below:

SALES CODE	DESCRIPTION	PROSTHETIC INTERFACE			Ø IM-PLANT (MM)	TRANSMU-COSAL HEIGHT (MM)	SCREW	
		MT	EH	IH			THREAD	CODE
TBM-0408	TI-BASE 4 CEREC® CM 0,8MM	X			3,5-4,0-5,0	0,8	M1,6	PF-1610
TBM-0415	TI-BASE 4 CEREC® CM 1,5MM	X			3,5-4,0-5,0	1,5	M1,6	PF-1610
TBM-0430	TI-BASE 4 CEREC® CM 3,0MM	X			3,5-4,0-5,0	3,0	M1,6	PF-1610
TBM-0445	TI-BASE 4 CEREC® CM 4,5MM	X			3,5-4,0-5,0	4,5	M1,6	PF-1610
TBE-3500	TI-BASE 4 CEREC® HE Ø3,5		X		3,5	0,6	M2,0	PPP-2085
TBE-3515	TI-BASE 4 CEREC® HE Ø3,5 X 1,5MM		X		3,5	1,5	M2,0	PPP-2085
TBE-3530	TI-BASE 4 CEREC® HE Ø3,5 X 3,0MM		X		3,5	3,0	M2,0	PPP-2085
TBE-3545	TI-BASE 4 CEREC® HE Ø3,5 X 4,5MM		X		3,5	4,5	M2,0	PPP-2085
TBE-4100	TI-BASE 4 CEREC® HE Ø4,1		X		4,1-5,0	0,6	M2,0	PPP-2085
TBE-4115	TI-BASE 4 CEREC® HE Ø4,1 X 1,5MM		X		4,1-5,0	1,5	M2,0	PPP-2085
TBE-4130	TI-BASE 4 CEREC® HE Ø4,1 X 3,0MM		X		4,1-5,0	3,0	M2,0	PPP-2085
TBE-4145	TI-BASE 4 CEREC® HE Ø4,1 X 4,5MM		X		4,1-5,0	4,5	M2,0	PPP-2085
TBI-3300	TI-BASE 4 CEREC® HI Ø3,3			X	3,3	0,6	M2,0	PPP-2085
TBI-3315	TI-BASE 4 CEREC® HI Ø3,3 X 1,5MM			X	3,3	1,5	M2,0	PPP-2085
TBI-3330	TI-BASE 4 CEREC® HI Ø3,3 X 3,0MM			X	3,3	3,0	M2,0	PPP-2085
TBI-3345	TI-BASE 4 CEREC® HI Ø3,3 X 4,5MM			X	3,3	4,5	M2,0	PPP-2085
TBI-4000	TI-BASE 4 CEREC® HI Ø4,0			X	4,0-5,0	0,6	M2,0	PPP-2085
TBI-4015	TI-BASE 4 CEREC® HI Ø4,0 X 1,5MM			X	4,0-5,0	1,5	M2,0	PPP-2085
TBI-4030	TI-BASE 4 CEREC® HI Ø4,0 X 3,0MM			X	4,0-5,0	3,0	M2,0	PPP-2085
TBI-4045	TI-BASE 4 CEREC® HI Ø4,0 X 4,5MM			X	4,0-5,0	4,5	M2,0	PPP-2085
SPM-4300	SCAN POST 4 CEREC® CM	X			3,5-4,0-5,0	5,2	M1,6	PF-1610
SPE-4300	SCAN POST 4 CEREC® HE		X		3,5-4,1-5,0	5,2	M2,0	PPP-2085
SPI-4300	SCAN POST 4 CEREC® HI			X	3,3-4,0-5,0	5,2	M2,0	PPP-2085

⚠ Note: For External Hexagon implants and Internal Hexagon diameter 5,0mm, the use of platform is indicated Switching.

**APPLICATION**

**TITANIUM BASE TIBASE 4 CEREC® EMFILS:** The Titanium Base TiBase 4 CEREC® Emfils is indicated according to the space Interocclusal Available, height of transmucosal and three-dimensional position of the implant.

They should be used in cementation-based unitary restorations on the titanium Base TiBase 4 CEREC® Emfils which are screwed on the implant. All interfaces are indicated for bolted unit restorations on dental implants. Check the compatibility between the chosen prosthetic interfaces. Titanium Bases are available in the Morse Taper (MT), External Hexagon (EH) and Internal Hexagon (IH) prosthetic interfaces, and are compatible with implants with the respective prosthetic interfaces. It is recommended to use the technique CAD/CAM technology through the Sirona Dental CAD/CAM System.

**SCANPOST 4 CEREC® EMFILS:** The Scanpost 4 CEREC® Emfils is indicated to facilitate scanning in hard to reach areas. The ScanPosts are available in the Morse Taper (MT), External Hexagon (EH) and Internal Hexagon (IH) prosthetic interfaces, and are compatible with implants with the respective prosthetic interfaces. It is recommended to use the technique CAD/CAM technology through the Sirona Dental CAD/CAM System.

**CONTRAINDICATIONS**

- The Titanium Base TiBase 4 CEREC® Emfils does not allow customization;
- It is contraindicated in restorations with angulation correction of more than 20 ° in relation to the implant axis and restorations whose length exceeds a ratio of 1:1.25 to Implant length (Example: Implant 10mm, not exceed 12,5mm);
- It is contraindicated for cases of multiple restorations;
- Not compatible with the Software Sirona Galileos®
- This product is contraindicated for patients presenting signs of allergy or hypersensitivity to the chemical components of the material: Titanium alloy Ti6Al4V-ELI (ASTM F136);

**MAGNETIC RESONANCE IMAGING (MRI) – INFORMATION**

These products are manufactured in metallic material that can be affected by MRI energy.

For more information, see “Security Information | Magnetic resonance imaging (MRI) “in <https://www.emfils.com.br/ifu>

**HANDLING**

FOR INTRAORAL SCANNING:

**TI BASE 4 CEREC® EMFILS:** Select and install the Titanium Base TiBase 4 CEREC® Emfils About the Implant according to previous planning. Enter the Transfer Scanning (Scanbody)\* Compatible With the equipment of the CAD/CAM system used on the titanium Base TiBase 4 CEREC® Emfils, and then run the Scanning according to manufacturer’s instructions Sirona. Make sure the correct fitting of the Scanbody\* About the Titanium Base TiBase 4 CEREC® Emfils.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select Tibase.

**SCANPOST 4 CEREC® EMFILS:** Select and install The Scanpost 4 CEREC® Emfils About the Implant according to previous planning. Enter the Transfer

Scanning (Scanbody)\* Compatible With the CAD/CAM system equipment used on the The Scanpost 4 CEREC® Emfils, and then run the Scanning according to manufacturer’s instructions Sirona. Make sure the correct fitting of the Scanbody\* About the Scanpost 4 CEREC® Emfils.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select ScanPost.

\* the Transfer Scanning (Scanbody) compatible with titanium Base Tibase 4 CEREC® Emfils is an accessory supplied by Dentsply Sirona.

**FOR GYPSUM MODEL SCANNING:**

TI BASE 4 CEREC® EMFILLS: Perform the molding of the implant according to the Previous planning, patient’s clinical situation and implant interface. Use artificial gums In the plaster model to simulate the soft tissue of the patient. Manufacture the plaster model according With the appropriate techniques. Select and install the Titanium Base TiBase 4 CEREC® Emfils On the implant analog. Enter The Transfer Scanning (Scanbody)\* Compatible with the equipment of the CAD/CAM system used and perform the scanning according to the manufacturer’s instructions. Make sure the correct Fitting the Transfer On the titanium Base TiBase 4 CEREC® Emfils.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select TiBase.

SCAN POST 4 CEREC® EMFILLS: Perform the molding of the implant according to the Previous planning, patient’s clinical situation and implant interface. Use artificial gums In the plaster model to simulate the soft tissue of the patient. Manufacture the plaster model according With the appropriate techniques. Select and install Scanpost 4 CEREC® Emfils on the implant analog. Enter the Transfer Scanning (Scanbody)\* Compatible with the equipment of the CAD/CAM system used and perform the scanning as shown in the of the manufacturer. Make sure the correct Fitting the Transfer About the Scanpost 4 CEREC® Emfils.

⚠ Note: After selecting the NBB 3.4 L Library, in the module of choice, select ScanPost.

\* the Transfer Scanning (Scanbody) compatible with titanium Base Tibase method 4 CEREC® Emfils is an accessory supplied by Dentsply Sirona.

**MESOSTRUCTURE DESIGN:**

Design the prosthesis structure in Sirona InLab Software (version 3.65) or Sirona CEREC® Software (version 4.2) According to previous prosthetic planning and the material to be used. Titanium Bases TiBase 4 CEREC® Emfils About implants Emfils (MT, EH and IH) Can be used in conjunction with the products Sirona, presented in the following table:

Note: the ScanPost Compatibility according to the table below:

⚠ Note: The Software does not permit the use of Scan Posts 4 CEREC® Emfils for jobs where Ti-Bases 4 CEREC® Emfils are used with 1,5mm / 3,0mm e 4,5mm transmucosal heights.

SIRONA PRODUCTS FOR USE IN CONJUNCTION WITH TIBASE 4 CEREC® EMFILLS				
LIBRARY	SCANBODY (SIRONA)	REF. SCANBODY OMNICAM (SIRONA)	REF. SCANBODY BLUECAM inEos X5 and inEos Blue (SIRONA)	CEREC BLOCKS RECOMMENDED
NBB 3.4L	L	6431329	6431303	inCoris ZI meso F0.5 (L) inCoris ZI meso F2 (L) IPS e.max CAD (L)

⚠ Note: The thickness of the minimum bone wall around the wall of the screw grooving of the structure is variable

PRODUCTS FOR USE IN CONJUNCTION WITH SCANPOST 4 CEREC® EMFILLS					
SCANPOST 4 CEREC® (EMFILLS)	LIBRARY	SCANBODY (SIRONA)	REF. SCANBODY OMNICAM (SIRONA)	REF. SCANBODY BLUECAM inEos X5 and inEos Blue (SIRONA)	TIBASE 4 CEREC® (EMFILLS)
SPM-4300	NBB 3.4L	L	6431329	6431303	TBM-0408
SPE-4300					TBE-3500 TBE-4100
SPI-4300					TBI-3300 TBI-4000

according to the material, as shown in the following table. The maximum angulation of the structure shall not exceed 20°. The conicity of the structure shall not exceed 6° (when Coping). In the case of an angled structure the height of the cembale portion from the emergency profile (prosthetic height), should not exceed 1:1.25 in ratio to Implant length.

MATERIAL	MINIMUM THICKNESS(MM)
ZIRCONIA	0,5
IPS e-max® CAD	0,9

- Be sure not to exceed a 20 ° angle between the implant shaft and the shaft of the restoration.
- If the Mesostructural is designed to receive the aesthetic ceramics, make sure that this will not narrow the screw grooving.
- The fitting cavity in the titanium Base Tibase 4 CEREC® Emfils mesostructural should not be coated.
- Make sure there are no sharp edges or corners in the design of the Mesostructural.
- Before making the structure of the prosthesis, make sure that the size of the machining block is compatible with the design to be carried out.
- Manufacture the workpiece as shown, using machining equipment Sirona CEREC® MC X, InLab MC XL or InLab MC X5.

**CLINICAL PROOF OF PROSTHESIS STRUCTURE:**

1. Perform a clinical proof of the structure Made on the titanium Base Tibase 4 CEREC® Emfils.
2. Make sure the proper fitting between The parts clinically and through radiography.
3. Make sure the occlusal space exists, Proximal Inter And that the structure of the prosthesis meets Aesthetic and functional requirements.
4. Finish the restoration according to the planning Appropriate techniques.

**CEMENTING OF THE TITANIUM BASE TIBASE 4 CEREC® EMFILLS:**

1. It is recommended to screw the titanium Base Tibase 4 CEREC® Emfils On an implant analog during The entire finish and finishing process of the prosthesis.
2. Position the Titanium Base Tibase 4 CEREC® Emfils On the implant analog and transcend the Fastening on the same with a slight tightening Using the key Installation indicated (see table below).
3. The surfaces of the titanium Base Tibase 4 CEREC® Emfils Aimed at cementing the structure of the Prosthesis should be treated with sandblasting (Aluminum oxide, 50 µm, maximum pressure of 2 bar). Protect screw access during The cementing process.
4. ⚠ Note: It is recommends the use of chemically resinous cements Activated for bonding in metals (e.g.: PANAVIA – Kuraray). In Disilicate Lithium (IPS E.max® CAD) is mandatory the use of IVOCLAR cement Multilink. The handling of the cement must follow the Manufacturer’s instructions.
5. Apply Cement to External portion of the titanium Base Tibase 4 CEREC® Emfils and press the reset by fitting according With the Indexer tab.
6. The restoration must be Pressed on the titanium Base Tibase 4 CEREC® Emfils And the excess of cement that is extravasated by the Orifice of the same should be immediately Removed.
7. After the cement’s prey Infrastructure is unscrewed

from implant analog and Remove excess cement remaining in the Titanium Base Tibase 4 CEREC® Emfils.

**CUSTOM PROsthESIS INSTALLATION:**

Sanitizer and sterilize the set, titanium Base Tibase 4 CEREC® Emfils and custom prosthesis, before To take it in the mouth. After these procedures, position the set in the mouth using the installation driver and the torque indicated as shown in the following table. Perform the occlusal adjustments and Interproximal needed.

SCREW MODEL	TORQUE (N.cm)	INSTALLATION DRIVER
PF-1610	25	HEXAGONAL 1,2MM
PPP-2085	30	HEXAGONAL 1,2MM

**PRESENTATION FORM**

All parts that make up the system 4 CEREC® Emfils are individually packaged in surgical paper packaging type Tyvek® and transparent film display. They are provided clean because they undergo a decontamination process, which eliminates any kind of dirt, however, require sterilization prior to clinical use.

**STERILIZATION**

The Tibase 4 CEREC® Emfils are supplied Non-sterile. All models are indicated for single use only.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization by Saturated Steam - Autoclaves

In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.

All surgical instruments and prosthetic components must be sterilized prior to use.

Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

**PRECAUTIONS**

- Do not use if The packaging is violated;
- Check the correspondence between The type of screw according to the implant model with its prosthetic interfaces;
- Before installing the products, make sure that they have the appropriate form of use for the selected prosthetic system;
- Special care must be taken during the project of the prosthesis so that the space Interocclusal No Insufficient or excessive.
- Ensure that the parts are not swallowed or aspirated by the person who is receiving the product (patient);
- All models of TiBases 4 CEREC® Emfils Are Indicated For single use only. Reuse of this product can Cause:
  - *Effects Biological adverse effects resulting from residual residues, be they, microorganisms and/or substances arising from previous uses;*

- *Changes in physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may impair the functionality of which it was intended;*
- *Reuse This product does not guarantee its safety and efficacy, and exempts any warranty from the related products.*

- Higher insertion torques than recommended or use of improper installation driver can cause damage to the material and make the system unusable;
- Always use the product sequence Emfils. The use of instrumental and/or Third-party components do not Ensures the perfect function of the sYstem of iMplante Emfils and exempts any warranty from the product.
- Improper planning can compromise the performance of the set Prosthesis, resulting in system failures such as loss or fracture of the implant or component, looting or fracture of the prosthetic screw;
- At the time of installation Of the component in the implant, Transgingival or implant analog, be sure to align The connection and/or o Screw, thus avoiding damage to the thread or locking of the system;
- The Titanium Base Tibase 4 CEREC® Emfils is not indicated for cementation in the mouth;
- Communicate to the laboratory that any modification or damage to the prosthetic interface requires that the workpiece be discarded and remade;
- When choosing titanium Base Tibase 4 CEREC® Emfils Ensure that it is compatible with the interface and diameter of the implant installed;
- The Titanium Base Tibase 4 CEREC® Emfils is not suitable for multiple restorations;
- Check the compatibility of Scanbodies, the With the Sirona software used;
- The components of the system 4 CEREC® Emfils Must Be perfectly seated to implant After applying the required torque to the Screw. This check is recommended using Periapical radiographs;

**ADVERSE EFFECTS**

Adverse effects will only occur if the choice and use of the product is inadequate.

**STORAGE CONDITIONS**

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

**DISPOSAL OF MATERIAL**

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

**DATE OF EXPIRATION**

Written on the product label.

CEREC® is a registered trademark of Sirona Dental Systems GmbH. Galileos is a mark used by the Sirona Dental Systems GmbH Germany. IPS E.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein. Multilink® is a registered trademark of Ivodiar Vivadent AG, Liechtenstein. Panavia™ is a registered trademark of Kuraray Co, LTD, JP.

The aforementioned libraries are protected by copyright and belong to their respective authors. All rights reserved to such authors. These libraries are available in Sirona InLab Software (version 3.65) or Sirona CEREC® Software (version 4.2) of the Sirona Dental Systems GmbH. All rights reserved. Not all products are available in all countries. Favor Contact Your authorized Distributor.

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