IFU

(Instruction for use)

Implante Corpo Único O'Ring / One Body O'Ring Implant / Implante Cuerpo Único O'Ring

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION

The Emfils one body o’ring implants are made of titanium (grade 5), according to ASTM F136. They are made with small external diameter because they are used temporarily. It has an external surface mechanically treated with a blasting operation with aluminum oxide, finished with acid attack. In its external macrogeometry, the one body abutment implant has an external thread that has variable pitch, which provides high primary stability. At the end of its cervical portion the implant has a ball-type connection that is designed for the mini-female (MFOR-4140) under pressure, in addition to a 2.0 mm external square geometry that serves as the interface for the installation driver.

The one body o’ring implant should be used in conjunction with the mini-female (MFOR-4140), spacer ring (EPRO-0001) and spacer capsule (PROC-0405) for temporary installation and retention of removable overdenture or partial dentures.

The nominal length of the implant is given by the distance from the apex of the implant to the groove below the square connection. The maximum installation torque for morse taper implants is 70 N.cm.

The Emfils one body o’ring implants are presented in the table below:

One Body Abutment Implant					
Sales Code	Diameter	Lenght	Sales Code	Diameter	Lenght
ICOE-2510	2,5	10,0	ICOE-3010	3,0	10,0
ICOE-2511	2,5	11,5	ICOE-3011	3,0	11,5
ICOE-2513	2,5	13,0	ICOE-3013	3,0	13,0
ICOE-2514	2,5	14,5	ICOE-3014	3,0	14,5
ICOE-2516	2,5	16,0	ICOE-3016	3,0	16,0

APPLICATION

The Emfils one body o’ring implant is indicated for surgical procedures in maxilla or mandible bones, providing temporary retention of removable overdenture or partial dentures.

In case of doubts about the indication, the use must be interrupted until all the questionable points have been clarified.

CONTRAINDICATIONS

It is not recommended to install the implant in the presence of acute infectious or inflammatory processes, inadequate bone volume or quality, serious medical problems such as disorders of bone metabolism, disorders of blood coagulation, inadequate healing, inadequate oral hygiene, incomplete jaw growth, patient uncooperative and unmotivated, drug or alcohol abuse, psychoses, prolonged functional disorders that resist any drug treatment, dry mouth, weakened immune system, diseases that require the regular use of steroids, endocrinological diseases, pregnancy or any other disorder, be it systemic or psychological, which may harm the tissue repair process.

MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION

The Emfils Implant System and components has not been evaluated for safety and compatibility

in the MR environment.

It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Emfils Implant System and components in the MR environment is unknown.

Scanning a patient who has this device may result in patient injury.

HANDLING

The dental implant should be used according all the requirements of biosafety and surgical asepsis.

The implant should be handled only in a sterile environment, by properly trained professionals and in appropriate attire at the time of surgery. The same must happen with the other instruments that will be used in the procedure.

All procedures must be performed by qualified professionals and the handling of the product is inherent in the training of the professional who will be responsible for using the components and instruments in the best possible way, respecting the Emfils indications for each product.

Drilling:

Once the surgical bed has been made, the implant must be removed from its packaging through the corresponding driver, taken to the

artificial alveolus and inserted by threading until the planned length.

Drill sequence recommended by Emfils in drilling for implant installation:

One Body Abutment Implant				
Drill Diameter	Bone Type			
	I and II		III and IV	
	Implant Ø		Implant Ø	
	2,5	3,0	2,5	3,0
2	X	X	X	X
2,5		X		X

In the bone bed drilling procedure, it is recommended to check the instruction for use of the Emfils surgical drills for more details on the application, handling and precautions of the instrument.

Implant installation:

Identify the planned implant model to be installed through the label on the back of the carton pack. The identification is given easily through text and colors, in the case of one body o’ring implants, the color is green. This label also identifies the diameter, length, batch number, date of manufacture and the validity of the implant. The diameter is also easily identified by the colors of the plastic tube cap on the front of the package: blue for implants of 2.5 mm diameter and yellow for implants of 3.0 mm diameter. At the top of this cap there is a circular label that identifies the length of the implant. In this front part of the package, we can also identify the circular label indicative of sterilization, which will appear in red color, if this label is in another color, or the product is not sterilized or has the sterilization expired and, in this case, the implant should not be used.

To remove the implant from the packing, remove the carton to access the internal blister and pull the tyvek surgical paper through the corner of the blister, accessing the inner tube containing the implant without the need of gloves. Lay the sterile tube over the surgical field without touching it.

With sterile surgical gloves, hold the bottle with a non-dominant hand and open the lid with the other hand, always keeping the tube upright, that is, with its opening facing up, preventing the implant from leaving the packaging.

In the opening of the cap, the user has access to the implant connection, which can be collected with the specific installation driver, searching the perfect fit between the driver and the implant by pressing lightly against each other. The installation driver for one body implants does not have a collection system that secures the implant to the driver, so after fitting the implant to the driver, the user must tuck the transparent tube through 90 ° so that the implant is removed from the implant packaging by gravity, without any manipulation other than that of the key itself. At this point the user must take care that the

implant does not fall, leaving it in an upright position for transport to the surgical socket. Transport the implant to the surgical alveolus and begin the installation by turning the driver clockwise. Complete the implant installation until planned final position. The installation driver must be connected to the Emfils torque wrench (TQ-0480), or a counter-angle motor. It must be ensured that the torque, during installation of the implant, does not exceed the recommended maximum torque of 70 N.cm.

PROSTHETIC SEQUENCE

In the immediate postoperative period, according to the chosen loading protocol, install the prosthesis respecting its indications and limitations.

PRE-OPERATIVE AND PROSTHETIC PLANNING

The model, diameter, length, position and quantity of implants should be selected for each clinical case, taking into consideration the anatomy, region to be rehabilitated, bone quality and quantity, and available space. Diagnostic closure of the clinical case should be performed, ensuring adequate alignment between implant (s), prosthetic component (s) and prosthesis (s). When installed in pairs, one body o’ring implants need to be installed perfectly parallel to each other, preferably bicortical.

PRESENTATION

The Emfils implants are conditioned in a package that offers triple protection: the primary package that contains the implant and is composed of titanium housing, plastic tube and flexible cap, the secondary package which contains the tube and it’s composed of blister sealed with tyvek surgical paper, and the tertiary package containing the sealed blister and is composed of folded carton paper and sealing by tamper proof label. This external label contains the identification and all information pertinent to the implant to be installed. Inside the tertiary carton, we also have the traceability label that is detachable on 3 smaller labels: one for the user to place on the patient’s medical record, another for the patient himself and the last label for reserve.

STERILIZATION

This product is indicated for single use and supplied sterile by the gamma radiation method, therefore, they are supplied ready to use.

PRECAUTIONS

- Do not use if the product packaging is broken;
- The implant sterilization is only guaranteed if the blister is not damaged;
- This product should be used immediately

after opening the blister, at the time of surgery. If you don’t use, dispose it;

- Do not use the implant with expired validity;
- Before installing the implant, make sure that it has the appropriate using way for the selected prosthetic system;
- Observe the intraoral tissues conditions, bone quality and bone quantity of the implant surgical bed through clinical, radiographic and / or tomographic exams. Failure to perform pre-surgical evaluation may compromise the success of the procedure;
- About the systemic aspect, consider the general health of the patient. In particular, care should be taken in cases of patients with drug allergies, local or systemic factors that may interfere with the healing processes of bone or soft tissues or in the process of osseointegration, for example: bone already exposed to radiation in the head or neck area, diabetes mellitus, anticoagulation drugs, hemorrhagic diathesis, bruxism, parafunctional habits, anatomically unfavorable bone situation, tobacco abuse, uncontrolled periodontitis, treatable maxillary pathologies and oral mucosa alterations;
- Bisphosphonate therapy results in potential risk of peri-implant osseonecrosis;
- The material to be used during the procedure must be sterile;
- Before each procedure, make sure that the parts fit together perfectly;
- The maximum installation torque is 70 N.cm. Insertion torque greater than that maximum, may can make the system inoperative;
- Ensure that the parts are not swallowed or aspirated by the person who is receiving the product (patient);
- This product is for single use and should not be reesterilized;
- The reusing of this product may cause:
 - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*
 - *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
 - *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as

loss or fracture of the implant, loosening or fracture of the mini-female (MFOR-4140);

- When installing components on the installed implant, be sure to align the connection with the mini-female (MFOR-4140) and always use the spacer ring (EPRO-0001), thus avoiding implant or mini-femoral damage;
- Before each procedure, check the conditions of the Emfils instruments, always respecting their useful life. Replace the instruments in case of damage, erasure, deformation and excessive wear;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Instruct the patient that if he or she notices that the prosthesis will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component;
- The clinical control of each prosthetic work executed with implants should be done periodically, so that the implant dentist can check the balance of the system or intervene clinically, if necessary;
- Although the implants or components are exceptionally resistant, they are subject to fractures in adverse biomechanical situations and / or inadequate rehabilitation plans;
- Transitional prostheses retained by one body o-ring implants are implant-retained-muco-supported and therefore should not transfer masticatory efforts to implants.

ADVERSE EFFECTS

The installation of implants, just like any other surgical procedure, can cause mild discomfort and localized edema. More persistent symptoms may occur as: chronic pain related to the dental implant, permanent paraesthesia, dysesthesia, loss of maxillary / mandibular bone border, localized or systemic infection, oro-antral or oro-nasal fistula, adjacent teeth adversely affected, irreversible damage to the adjacent teeth, fractures on implant, maxillary, bone or prosthesis, prosthetic problems, nerve damage, exfoliation, hyperplasia.

The eventual failure of the osseointegration process, or even after osseointegration, with the loss of the implant and the prosthesis can occur due to: inadequate osteotomy, infections, diseases or systemic problems, decrease in remaining bone quantity or quality, parafunction, lack or failure of irrigation, use of non-specific instruments and / or no cutting power, poor oral hygiene, lack of cooperation on the part of the patient , occlusal trauma, lack of prosthetic passivity and absence of specific training.

PÓS-OPERATIVE CARE AND MAINTENANCE

Orient the patient about the need to perform

a professional follow-up after surgery and obey the care guidelines, hygiene of the surgical area and prosthesis, as well as the prescription of medications. These guidelines are the responsibility of the professional.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

All products and consumables used in a surgery for the implant dental installation can endanger the health of those who handle them after use. Before disposing them in the environment, it is recommended to consult and comply with the current regulations.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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