IFU

(Instruction for use)

Transferente e Coping Plástico / Impression and Plastic Coping / Transfer y Casquillo Plástico

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfils instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

DESCRIPTION
Transfers:

Transfers are part of the prosthetic components family of transitory use to transfer the exact three-dimensional position of the implant or base abutment that is installed in the mouth. This transfer can be performed through two techniques: open tray or closed tray, so the Emfils transfers enable the use of the technique that best meets the needs of the professional, as they contemplate both techniques. The product is available with different fittings, depending on the type of implant or base abutment used.

The raw material used in the transfers is stainless steel surgical series 300, according to ASTM F899 303.

Plastic Copings:

The plastic copings are part of the prosthetic components family of the Morse taper line of transitory use or even definitive, if they are waxed and casted. The copings are only indicated if you are using the abutments for cementing of the Morse taper line. They can be used to transfer the exact three-dimensional position of the straight or angled abutment that is installed in the mouth or used in the laboratory operation of waxing and casting, forming the structure of the prosthesis that is cemented in the straight or angled abutment, because the interface with these products will be maintained. The copings are available with the exact fit of the straight and angled pillars of the Morse taper system and are made of calcifiable polymer.

The transfers are presented in the table below:

Código de Venda	Tipo Conexão								Moldeira		Interface		Diâmetro Implante ou Transmuc.	Parafuso
	CM	LZ	HE	HI	Transmucoso				Aberta	Fechada	Rotac.	Não Rotac.		
					Multi	Mini 4,0-5,0	Mini 4,8	Angul.						
CPTMA-3312				X					X			X	3,3-3,5	M2
CPTMA-4012				X					X			X	4,0	M2
CPTMA-5012				X					X			X	5,0	M2
CPTMF-3312				X						X		X	3,3-3,5	M2
CPTMF-4012				X						X		X	4,0	M2
CPTMF-5012				X						X		X	5,0	M2
CTAE-3512			X						X			X	3,5	M2
CTAE-4012			X						X			X	4,0	M2
CTAE-5012			X						X			X	5,0	M2
CTAI-3512		X							X			X	3,5	M1.6
CTAI-4012		X							X			X	4,0	M1.6
CTAI-5012		X							X			X	5,0	M1.6
CTAM-3512	X								X			X	3,5-4,0-5,0	M1.6
CTFE-3512			X							X		X	3,5	M2
CTFE-4012			X							X		X	4,0	M2
CTFE-5012			X							X		X	5,0	M2

Código de Venda	Tipo Conexão								Moldeira		Interface		Diâmetro Implante ou Transmuc.	Parafuso
	CM	LZ	HE	HI	Transmucoso				Aberta	Fechada	Rotac.	Não Rotac.		
					Multi	Mini 4,0-5,0	Mini 4,8	Angul.						
CTFI-3512		X								X		X	3,5	M1.6
CTFI-4012		X								X		X	4,0	M1.6
CTFI-5012		X								X		X	5,0	M1.6
CTFM-3512	X									X		X	3,5-4,0-5,0	M1.6
CTNAT-4012					X				X			X	4,0	M2
CTNAT-5012					X				X			X	5,0	M2
CTNFT-4012					X					X		X	4,0	M2
CTNFT-5012					X					X		X	5,0	M2
CTRAM-4012						X			X		X		4,0	M1.6
CTRAM-5012						X			X		X		5,0	M1.6
CTRAT-4012					X				X		X		4,0	M2
CTRAT-5012					X				X		X		5,0	M2
TRFM-4012						X				X	X		4,0	M1.6
TRFM-5012						X				X	X		5,0	M1.6
TRFT-4012					X					X	X		4,0	M2
TRFT-5012					X					X	X		5,0	M2
TFMP-4810							X	X		X	X		4,8	M1.4
TMPA-4812							X	X	X		X		4,8	M1.4

The copings are presented in the table below:

Código de Venda	Design		Comprimento		Pilares de Utilização
	Regular	Estético	Curto	Longo	
CCE-4506		X	X		CECM-3515 - CECM-3530 - CECM-3545
CCR-3506	X		X		CRCM-3515 - CRCM-3530 - CRCM-3545 CACM-3515 - CACM-3530 - CACM-3545
CLE-4508		X		X	CELM-3515 - CELM-3530 - CELM-3545
CLR-3508	X			X	CRLM-3515 - CRLM-3530 - CRLM-3545 CALM-3515 - CALM-3530 - CALM-3545


APPLICATION

Transfers:


The transfer allows to transfer, through molding, the three-dimensional position of the implant or base abutment installed in the mouth. The solution is presented for closed and open tray molding techniques. The transfers are of transitory use and must be installed in the mouth, on the implant or base abutment, only during the molding process. After the molding is done, the transfers must be removed and used in the plaster model manufacturing process.

The transfer has a geometry that allows the transfer of the three-dimensional position of the implant or base abutment through the technique of the molding. The product may have a rotational or non-rotational fit, depending on the interface and the type of prosthesis to be used.

To fix the transfer, you must use the corresponding screw with the system and the technique that is being used.

Transfer to closed tray technique:

Copings:

The copings allow to transfer, through molding, the three-dimensional position of the implant or base abutment installed in the mouth. The copings are of transitory use, but may become definitive if they are used in the laboratory operation of waxing and casting, forming the structure of the prosthesis that is cemented in the straight or angled abutment of the Morse taper system, because the interface with these products will be maintained.

CONTRAINDICATIONS

Not suitable for multiple transfer of several implants where there is not good parallelism between them, because it may occur if there is no withdrawal axis.

MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION

The Emfils Implant System and components has not been evaluated for safety and compatibility in the MR environment.

It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Emfils Implant System and components in the MR environment is unknown.

Scanning a patient who has this device may result in patient injury.

HANDLING

Transfer to open tray technique:


ATTENTION: The figures are merely illustrative. They do not represent the actual dimensions and characteristics of the product.

Transfer to open tray technique:

- a) Fit the transfer body in the implant or base abutment and fix it with the corresponding screw, and only turn the screw with the hexagonal driver that best suits the work. Do not exceed the recommended maximum torque.
- b) For multiple restorations, the union of the transfers with acrylic resin must be done;
- c) The tray must have an access opening for the transfer and the screw, this opening must be made if it does not have;
- d) Handle the molding material according to the manufacturer's recommendations;
- e) Load the tray with the material and take it to the patient's mouth;
- f) After allowed to set off the material, with the help of the hexagonal driver, loosen the screw;
- g) Remove the tray and make sure that the transfer has accompanied the mold;
- h) Carefully fit and screw the implant analogue onto the transfer;
- i) Make sure not to move the transfer during the fitting of the analogue;
- j) Make the plaster model according to the appropriate techniques;
- k) Remove the transfer of the plaster model.

Transfer to closed tray technique:

- a) Fit the transfer body in the implant or base abutment and fix it with the corresponding screw, and only turn the screw with the hexagonal driver that best suits the work. Do not exceed the recommended maximum torque.
- b) Make sure that the transfer is well seated;
- c) Handle the molding material according to the manufacturer's recommendations;
- d) Load the tray with the material and take it to the patient's mouth;
- e) After allowed to set the material, remove the tray;
- f) Using a hexagonal driver, loosen the screw and remove the transfer;
- g) Then screw the transfer into the analog of the implant or base abutment;
- h) Carefully fit the transfer-analog assembly into the mold so that a precise fit is obtained in the same position as in the mouth;
- i) Make the plaster model according to the appropriate techniques;
- j) Remove the transfer of the plaster model.

Transfer with coping:

- i) Fit the body of the coping by pressing it to the abutment so that it is well seated;
- j) Handle the molding material according to the manufacturer's recommendations;
- k) Load the tray with the material and take it to the patient's mouth;
- l) After allowed to set the material, remove the tray;
- m) Carefully insert the analog of the abutment into the coping so that it is well seated;
- n) Make sure not to move the coping during the insertion of the abutment analog;
- o) Make the plaster model according to the appropriate techniques;
- p) Remove the coping of the plaster model.

PRESENTATION

The transfer and copings are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

PRECAUTIONS

- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Ensure that instruments are not swallowed or aspirated by the person receiving the product (patient);
- Transfers and copings are for single use only, the reuse of this product does not

guarantee its safety and effectiveness and disclaims any warranty of the related products;

- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of Emfils products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfils from offering any product warranty;
- Special care should be taken during prosthesis design for what the interocclusal space is not insufficient or excessive;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;
- The component must be perfectly seated to the implant, transmucosal or abutment after application of the necessary torque on the screw or manual pressing;
- Open tray technique is not recommended for patients with poor oral opening;
- Not suitable for multiple transfer of several implants where there isn't good parallelism between them, because it may occur if there is no withdrawal axis.
- Too much attention to repositioning the transfer in the mold in the closed tray technique;
- In the closed tray technique, obliterate the screw hexagon with utility wax so that it isn't copied into the molding and makes it difficult to reposition the transfer in the mold.

ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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