## IFU

(Instruction for use)

Transmucoso Mini Angulado 4.8 mm / Angled Mini Base Abutment 4.8 mm / Pilar Cónico Mini Angulado 4.8 mm

PORTUGUÊS

ENGLISH

ESPAÑOL

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Favor contactar seu distribuidor autorizado.

This product has been designed to integrate clinical, mechanical and biological concepts, so, for best results it should only be used with original associated components according to Emfil instructions and recommendations.

All procedures must be performed by qualified professionals such as Implantodontists, Periodontists, Prosthodontists and Oral Maxillofacial Surgeons.

## DESCRIPTION

The angled mini base abutments are part of the prosthetic components family and are used as an intermediate between the implant and the prosthesis (crown). The product is available with different fittings, depending on the type of implant used, also has different collars that correspond to the gingival heights, taking into account the mucosal thickness variations. In its upper part it presents a conical geometry without hexagon that allows the connect of transfer, provisional, protector and rotational burnout abutment indicated for multiple prosthesis, through fixation by prosthetic screw.

The raw material used is the titanium alloy (grade 5), according to ASTM F136.

The angled mini base abutments are presented in the table below:

Sales Code	Connection Type				Piece		Implant Diameter	Angulation	Collar Height	Screw
	MT	ZL	EH	IH	Single	Compound				
CTAPE-31720			X			X	3,5	17°	2,0	M2
CTAPE-31740			X			X	3,5	17°	4,0	M2
CTAPE-33020			X			X	3,5	30°	2,0	M2
CTAPE-33040			X			X	3,5	30°	4,0	M2
CTAPE-41720			X			X	4,0	17°	2,0	M2
CTAPE-41740			X			X	4,0	17°	4,0	M2
CTAPE-43020			X			X	4,0	30°	2,0	M2
CTAPE-43040			X			X	4,0	30°	4,0	M2
CTAPI-31720				X		X	3,3-3,5	17°	2,0	M2
CTAPI-31740				X		X	3,3-3,5	17°	4,0	M2
CTAPI-33020				X		X	3,3-3,5	30°	2,0	M2
CTAPI-33040				X		X	3,3-3,5	30°	4,0	M2
CTAPI-41720				X		X	4,0	17°	2,0	M2
CTAPI-41740				X		X	4,0	17°	4,0	M2
CTAPI-43020				X		X	4,0	30°	2,0	M2
CTAPI-43040				X		X	4,0	30°	4,0	M2
CTAPM-1715	X					X	3,5-4,0-5,0	17°	1,5	M1.6
CTAPM-1730	X					X	3,5-4,0-5,0	17°	3,0	M1.6
CTAPM-3015	X					X	3,5-4,0-5,0	30°	1,5	M1.6
CTAPM-3030	X					X	3,5-4,0-5,0	30°	3,0	M1.6

## APPLICATION

The angled mini base abutments are indicated for cemented unit prostheses on the implant, according to available interocclusal space, existing gingival height and three-dimensional position of the implant. It can be used in immediate or late rehabilitation procedures, in the maxilla or mandible. Before finishing the clinical work, make sure that the screw has been tightened with the necessary torque.

## CONTRAINDICATIONS

The angled mini base abutments are not indicated for personalization. Contraindicated in cases where the interocclusal space is insufficient and / or the three-dimensional position of the implant is unsatisfactory.

## MAGNETIC RESONANCE IMAGING (MRI) – SAFETY INFORMATION

The Emfil Implant System and components has not been evaluated for safety and compatibility in the MR environment.

It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Emfil Implant System and components in the MR environment is unknown.

Scanning a patient who has this device may result in patient injury.

## HANDLING

After choosing the best abutment collar that can meet the client's gingival thickness, install the angled mini base abutment over the

implant, respecting the correct connection and applying the recommended torque. Always protect the access thread at the top of the part, preventing undesirable substances from obstructing the thread of the prosthetic screw. On the angled mini base abutment it's possible the fit of transfer, provisional, protector or burnout abutment with the same connection, fixing them by prosthetic screw.

## PRESENTATION

The angled mini base abutments are indicated for single use and individually packaged in tyvek® surgical paper and transparent film display. They are supplied cleaned because they undergo a decontamination process that eliminates any kind of dirtiness, but need to be sterilized before using.

## STERILIZATION

Supplied non-sterile. Indicated for single use.

STERILIZATION: is the procedure that aims at the total elimination of microorganisms, whether in vegetative or sporulated form.

Sterilization By Saturated Steam - Autoclaves

- In a conventional autoclave, the instrument should stand for 30 minutes at a temperature of 121 ° C when this temperature is reached.
- All surgical instruments and prosthetic components must be sterilized prior to use.
- Parts that have rubbers and plastics, for connection to another part, should not be autoclaved or sterilized being connected to each other. Example: Drill connected to Protractor.

## PRECAUTIONS

- Before installing the product, make sure it has the same prosthetic interface as the implant;
- Ensure that the stability of the implant is sufficient to support the installation torque of the prosthetic component and functional load, see instructions for use of the implant.
- Do not use if the product packaging is broken;
- Check the correspondence between the threads and prosthetic interfaces of the components;
- Ensure that the pieces are not swallowed or aspirated by the person receiving the product (patient);
- Any type of base abutment is for single use only;
- Reusing this product may cause:
  - *Adverse biological effects from remaining residues, whether microorganisms and / or substances from previous uses;*

- *Changes in the physical, mechanical and chemical characteristics of macro or micro structural order, altering the original conditions of the product that may harm from its intended functionality;*
- *The reuse of this product does not guarantee its safety and efficacy, and disclaims any warranty of related products.*
- Be sure to use a prosthetic screw compatible with the prosthetic interface and the prosthetic component;
- Insertion torques higher than recommended or use of improper connection drivers can cause damage to the material and make the system unusable;
- Always use the sequence of Emfil products. The use of instruments and / or components from other manufacturers and systems does not ensure the perfect function and exempts Emfil from offering any product warranty;
- Special care should be taken during prosthesis design for what the interocclusal space is not insufficient or excessive;
- Inadequate planning may compromise the performance of the implant / prosthesis assembly, resulting in system failure, such as loss or fracture of the implant or component, thread loosening or fracture of the prosthetic screw;
- At the time of installation, be sure to align the screw along the component shaft, thus avoiding damage inside the implant or locking the system;
- The component must be perfectly seated on the implant after the necessary torque has been applied to the screw. This check is recommended using periapical radiographs with the parallelism technique, when applicable.
- The angled mini base abutment is not recommended for use in the posterior regions of the mouth when the implant has a narrow diameter;
- The angled mini base abutment should be used as recommended, do not create additional angulation during customization;
- The selection of the material of the prosthesis structure should take into account local aspects of the patient;
- It is recommended that the final prosthesis be sanitized before installation in the mouth;
- Instruct the patient that if he or she notices that the healing cup will loosen or even detach from the position of the surgery due to out-of-control conditions, he or she should immediately go to the clinic to reattach or replace the component.

#### ADVERSE EFFECTS

Adverse effects will only occur if the choice and use of the product is inadequate.

#### STORAGE CONDITIONS

The product should be stored in a dry place and in the original packaging at room temperature and not expose to direct sunlight. Incorrect storage may influence the original characteristics of the packaging and the product.

#### DISPOSAL OF MATERIAL

The disposal of the material must comply with local regulations and environmental requirements, taking into account the different levels of contamination and selective collection.

#### DATE OF EXPIRATION

Written on the product label.

Not all products are available in all countries.

Please, contact the authorized distributor.

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